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Photo Guide for

LudgerClean™ S  
Glycan Cleanup Cartridges

Product # LC-S-A6

**Ludger Document # LC-S-A6-Photo-Guide-v1.0**

**Ludger Ltd**

Culham Science Centre  
Oxford OX14 3EB  
United Kingdom

Tel: +44 1865 408 554  
Fax: +44 870 163 4620  
Email: [info@ludger.com](mailto:info@ludger.com)  
[www.ludger.com](http://www.ludger.com)

## Photo Guide for LudgerClean™ S Glycan Clean Up Cartridges

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This is a step-by-step visual guide to using our LudgerClean S glycan clean up cartridges (Ludger Cat# LC-S-A6). You would typically use these for purification of labeled glycans after fluorescent tagging e.g. 2-AB, 2-AA or procainamide. This document complements the full [LC-S-A6 product guide](#).

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The sample to be cleaned must have a volume of 15 µl or less. If your sample has a greater volume then dry it down using a centrifugal evaporator and reconstitute in not more than 15 µl of water

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
Prime the LudgerClean™ S cartridges (one per sample) by washing each with the following :

- 1st wash - 1 ml water
- 2nd wash - 5 ml 30 % acetic acid (aq)
- 3rd wash - 1 ml acetonitrile

*Allow each wash to drain completely before adding the next. If flow is restricted, e.g. by an air gap, apply a slight pressure to the top of the cartridge in order to resume normal flow.*

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
Spot each sample onto a freshly washed cartridge disc ensuring that the disc is still wet with acetonitrile. Spread the spot over the entire disc surface if possible as this aids cleanup.

### 4 让样品吸附到薄膜上 Allow sample to adsorb onto membrane

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Allow adsorption for 15 minutes.


### 5 从样品瓶中添加剩余样品 Add residual sample from sample vial

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Rinse each sample vial with 100  $\mu$ l acetonitrile and apply to the corresponding cartridge disc.

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Wash each disc with 1 ml acetonitrile, followed by 5 x 1 ml 96 % acetonitrile\* / 4% water

*\*(for cleanup of O-glycans or N- and O-glycans labeled with procainamide, substitute with 100% acetonitrile)*

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Elute glycans off membrane into a suitable container by eluting with 2 x 0.5 ml water.

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Allow each 0.5 ml aliquot to drain before the next is applied.

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If appropriate, evaporate the glycan containing fraction to dryness, then dissolve in a desired volume of water or solvent for further analysis.

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Your glycans are now ready to analyse.